

Key Tips when using Tappy, Tok-Tok

1. It is highly recommended to **watch the Manual Video** before use.
(Video #48 can be **DOWNLOADED** at www.tappy.co.kr)
(Download at the portal site by searching [Tappy Tok-Tok])

2. Follow the procedures in the enclosed manual exactly.

Procedure: **1. FGF & EGF 2. nano Platinum Serum 3. Mask Pack**
Or: **1. EGF Serum 2. Platinum Serum 3. Mask Pack**

3. **Pre-tapping the skin 30-50 times maximizes stimulation and regeneration of skin.**

The first 30-50 times of tapping around eyes and the whole face is the most important step for stimulation of cell regeneration in the dermis.

During this step, the active constituent of serum (solution) can be delivered into the dermis through the **screw thread** on **each micro needle slowly**.

4. For safe use, **do not share** with other people.

The needle of Tappy is susceptible to contamination if it is stained with blood in the dermis. In order to protect from possible contamination only use **one Tappy per person**.

5. Perform this Therapy **2 times a week**. (8 times/month)

Tapping with Tappy punctures and forms very fine micro holes on the skin.

Through these holes, the active ingredient of serum is transferred into the dermis and the skin cells around the holes regenerate to granulate.

3-4 days is needed for Cell Regeneration, applying twice a week is the most effective way to achieve results.


6. **Apply Tappy on the skin only.**

The needle is highly durable, tested 30,000 times during tapping tests.

The micro needle is as thick as hair and has a subtle screw thread design which can be easily broken by dropping or tapping the needle on places other than skin.

Be careful to always store this device with its cap on.

7. **Please do not spin too hard** when combining bottle with Tappy(Tok-Tok) head, which might lead to cracks on them. You can use the product with **Tappy head lightly locked**.



If you have any questions using Tappy, feel free to contact belows
U-BioMed Inc. 409, BI Center, Kyungpook National University, Buk-Gu, Daegu, 702-701, KOREA

Mail: ubiomed@hanmail.net Tel.: +82-53-942-1117

H.P.: www.ubiomed.co.kr Mall: www.tappy.co.kr